## HYDROFERA BLUE® BACTERIOSTATIC DRESSINGS

### INSTRUCTIONS FOR USE

#### DESCRIPTION

Hydrofera Blue® Bacteriostatic Dressing is a sterile absorptive foam dressing made of bacteriostatic polyvinyl alcohol (PVA) sponge, Methylene Blue (less than or equal to 0.00025 g/g), and Crystal Violet (less than or equal to 0.00025 g/g).

#### INDICATIONS

Hydrofera Blue® Bacteriostatic Dressings are intended as an external dressing for use in local wound management, such as pressure ulcers, venous stasis ulcers, arterial ulcers, donor sites, abrasions, lacerations, superficial burns, postsurgical incisions, diabetic ulcers, and other wounds caused by trauma.

Hydrofera Blue® Dressings are appropriate for use with enzymatic debriders.

#### CONTRAINDICATIONS

Hydrofera Blue® Dressings are not indicated for third-degree burns.

#### APPLICATION OF HYDROFERA BLUE® DRESSINGS

##### 1) Wound Site Preparation

- Prior to the initial application of Hydrofera Blue® Dressing, the wound should be debrided of any necrotic tissue.
- The wound should be cleansed with normal saline and/or an appropriate cleansing solution.

##### 2) Dressing Application

- Open the Hydrofera Blue® Dressing package, keeping the dressing on the inside of the package.
- If the wound presents with surrounding erythema or rolled wound edges (epibole), the dressing should be larger than the wound to cover the periwound area.
- In other wounds, the dressing may be cut to fit the wound bed. An appropriate skin prep may be used in the periwound area.
- MOISTEN THE DRESSING WITH STERILE SALINE OR STERILE WATER, AND SQUEEZE OUT THE EXCESS.
- Apply the dressing to the wound. The dressing may be folded or layered to ensure desired fit.
- Secure with a secondary dressing to prevent displacement and to maintain moisture content.

Examples of secondary dressings include:

Level of Exudate	Secondary Dressing
Low	Semi-occlusive film Occlusive film Vaseline gauze
High	Bordered gauze Rolled gauze

#### Instructions for Dressings with Moisture Retentive Film

Prepare the wound site as described in step 1 above.

Dressing Application:

- Remove the dressing from the package, and hydrate the exposed foam side with approximately 25-30 cc of sterile saline or sterile water. Squeeze out excess sterile saline or sterile water.
- Place the hydrated foam side directly to the wound with the **film side on top**.
- Secure with a secondary dressing (see instructions in step 2f above).

#### Instructions for Island Dressings

Prepare the wound site as described in step 1 above.

Dressing Application:

- Remove the dressing from the package, and peel off the dressing backing.
- Apply the dry dressing to the wound, and lightly press adhesive to the skin.
- Using a syringe, hydrate the dressing with sterile saline or sterile water through the portal opening on top of the dressing.
- Remove the outer frame backing from the dressing, and ensure the adhesive is secured to the skin.

#### Instructions for Tunneling Dressings

Prepare the wound site as described in step 1 above.

Dressing Application:

- Remove the tunneling dressing from the package, and moisten the tip with sterile saline or sterile water to aid with ease of insertion. It is important to keep the packing rigid for placement.
- Gently insert the moistened end of the tube (softened) into the wound canal. If necessary, a rigid swab may be used to help guide it into place.
- After insertion, add additional sterile saline or sterile water to the end of the dressing so the entire tube is moist.
- Secure with a secondary dressing (see instructions in step 2f above).

#### Instructions for Ostomy Dressings

Prepare the wound site as described in step 1 above.

Dressing Application:

- Trim the inner diameter of Hydrofera Blue® Dressing to fit around the stoma site (as needed).
- Cut the outer diameter of the dressing slightly smaller than the ostomy pouch adhesive area (as needed).
- The ostomy adhesive should adhere directly to the patient's skin.

#### DRESSING CHANGE AND REMOVAL

1) The first dressing change should occur **at 24 hours**.

2) Examine the area of the dressing in contact with the wound bed:

- If the dressing has turned white or lightened in color, the dressing should be changed again at 24 hours until the dressing retains its color.
- If the dressing has retained its color where it is in contact with the wound, a new dressing can be applied and left in place for up to 72 hours.

3) **Do not allow the dressing to completely dry out.** The dressing should be rehydrated as needed. If a Hydrofera Blue® Dressing is allowed to dry out, rehydration with sterile saline or sterile water is recommended prior to removal.

4) The dressing should be changed if it turns white or upon strikethrough drainage.

#### PRECAUTIONS AND OBSERVATIONS

1) This material is not an implant.

2) **Infections:** All wounds can be expected to contain some microorganisms; however, if true clinical symptoms should develop, appropriate steps as defined by the attending physician should be taken to address that infection. Regular evaluation and cleansing of infected wounds should be a common practice.

3) **Underlying Causes:** Some wounds, such as leg ulcers and pressure ulcers, can be nonhealing unless steps are taken to correct the underlying pathology. Hydrofera Blue® Dressings are designed for local wound management and to be part of the overall management program for these types of wounds.

4) **Drying Out:** Hydrofera Blue® Dressings are designed to maintain a moist wound environment that will facilitate wound healing. If a Hydrofera Blue® Dressing is allowed to dry out, rehydration with sterile saline or sterile water is recommended prior to removal.

**Note:** Hydrofera Blue® Dressings have two organic pigments bound to the foam to render it bacteriostatic. The dressing has been tested *in vitro*, and found to be a bacteriostatic barrier for the following bacteria and yeast: Methicillin Resistant *Staphylococcus aureus* (MRSA), *Staphylococcus aureus*, *Staphylococcus epidermidis* (coagulase neg.), *Pseudomonas aeruginosa*, *Escherichia coli*, *Enterococcus faecalis*, Vancomycin Resistant *Enterococcus* (VRE), *Bacillus subtilis*, *Yersinia enterocolitica*, *Serratia marcescens*, *Candida albicans*.

The dressing has not been subjected to a double-blind clinical study.

All dimensions shown represent the product in its hydrated state.

All safety and effectiveness testing is on file at Hydrofera, LLC.

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Fabriqué par Hydrofera, LLC.


## HYDROFERA BLUE® DRESSING APPLICATION

**Remove** Hydrofera Blue® Dressing from the package, and thoroughly **moisten** it with **sterile saline** or **sterile water**.

Once Hydrofera Blue® Dressing is thoroughly moistened, **squeeze** out the excess liquid.

**Apply** the dressing to the wound. If the wound presents with the surrounding erythema or rolled wound edges, the dressing should be **larger** than the wound to cover the **periwound** area. In other wounds, the dressing may be **cut to fit** the wound bed.


**Cover** Hydrofera Blue® Dressing with a **secondary** dressing to prevent displacement and to maintain **moisture** content.

The **first** dressing change should occur at **24 hours**. Examine the area of the dressing in contact with the wound bed.

- If the dressing has retained its blue color where it is in contact with the wound, apply a new dressing and leave in place for up to 72 hours.
- If the dressing has turned from blue to white or lightened in color, the dressing should be changed every 24 hours until it retains its blue color. The dressing should be changed if it turns white or upon strikethrough of drainage.

Do not allow the dressing to completely **dry out**. The dressing should be rehydrated as needed. If the Hydrofera Blue® Dressing dries out, rehydration with sterile saline or sterile water is recommended prior to removal.

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## PANSEMENT HYDROFERA BLUE® APPLICATION

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Tous les tests d'innocuité et d'efficacité sont disponibles auprès de Hydrofera, LLC.